

MR Uyumlu Ses ve Müzik Sistemi

MRI Compatible Sound System

- Her marka ve model MR cihazıyla birlikte kullanılabilir.
- MR ihtiyaçlarına uygun, özel amfi sistemine sahiptir.
- Dahili amfi hoparlörleri sayesinde hem teknik oda, hem MR odasına müzik yayını yapılabilen, ses seviyeleri ayrı ayrı ayarlanabilmektedir.
- Hastayla iletişim kurabilmek için bas-konuş mikrofoni mevcuttur. MUTE devresi sayesinde, mikrofona konuşma anında müzik yayını otomatik olarak kesilir.
- fMRI uygulamalarında kullanım için uygundur.
- Bilgisayar, DVD çalar, MP3 çalar, USB veya herhangi bir kaynaktan müzik yayını yapılabilir.
- Ses, kafa koili içerisine sığabilecek incelekte, akustik kulaklıkla hastaya iletilmektedir.
- 3 Tesla'ya kadar MR uyumluluk sertifikasına sahiptir.

- Can be used with all MRI scanners.
- Specifically designed sound control unit for MRI needs.
- Integrated speakers in music control unit provide the same sound with patient in control room.
- Sound levels in MRI room and control room can be adjusted separately.
- Push-to-talk microphone is provided to communicate with patient.
- MUTE circuit pauses music when the microphone is active.
- System is suitable with fMRI applications.
- PC, DVD player, MP3 player, USB or any other source can be used.
- MRI compatibility (up to 3 Tesla) certificate.

